


# SNOREKIL

## HOW IT WORKS

### SNOREKIL MAKES YOU STOP SNORING BY

gently repositioning your lower jaw in order to open the airway at the back of the throat (fig. 1)

Pushing your lower jaw forward (fig. 2) gently stretches the soft tissue at the back of the throat



SNOREKIL comes in two parts (fig. 1) Top mouthguard (fig. 2) The bottom mouthguard (2) has the "gate" which allows you to adjust the effect of the device on your lower jaw position



Tel: 01492 541886

# SNOREKIL

has been developed by our in-house dental professional with over 35 years experience.

## WHY IS SNOREKIL BETTER THAN OTHER STOP SNORING DEVICES

### Snorekil

Easy to use, 2 simple parts, just click together  
Fully adjustable - simple gate system, soft plastic fits to shape of your mouth.  
Low cost - affordable oral device, less than £50, no need for special fitting  
Comfortable to wear at night - breath through mouth and nose - suitable for hay fever sufferers

### Other stop snoring devices

Can be complicated to put together  
May be hard to fit and adjust, may need special fitting by a dental professional  
Can be expensive - some more than double the price of Snorekil  
Some oral devices do not allow you to breath through your mouth

For more information visit

[www.snorekil.co.uk](http://www.snorekil.co.uk)  
[www.snorekil.com](http://www.snorekil.com)

SleepWorks Products Ltd  
presents

# SNOREKIL

## Anti-Snoring Dental Device

AWARD WINNER AT THE  
INTERNATIONAL EXHIBITION OF INVENTORS GENEVA 2007  
FOR ITS INNOVATIVE DESIGN

WINNER AT THE  
BRITISH  
INVENTION



OF THE YEAR  
AWARDS

AWARD WINNER  
AT



MEDI WA  
2007

Tel: 01492 541886